

AT-1 Applanation Tonometer

Model name	AT-1 Applanation Tonometer	
Attachable slit lamp models	700GL, 700GL NSW, 700GL LPH, 4ZL, 2ZL	30GL
Mounting tonometer base	TB700	TB30
Measurement range	0-80 mmHg	
Minimum scale value	2 mmHg	
Measurement method	Goldmann method (Goldmann Applanation Tonometry)	
Measurement deviation	Nominal value $\pm 1.5\%$ or ± 0.49 mN, whichever greater	
Reverse span (Hysteresis)	≤ 0.29 mN	
Diameter of applanation field	3.06 mm	
Depression area	7.354 mm ²	
Dimensions	72 (W) x 135 (D) x 251 (H) mm	
Weight	0.44 kg	

*Specifications as detailed in this sheet may be subject to change as improvements are made to products.

*For the latest information, please visit our website (<https://www.takagi-j.com/en/>).