


# AT-1 / AT-2

IOP measurement for  
routine clinical practice

- Accurate and stable measurements based on the Goldmann applanation tonometry method
- Seamless prism design for easy cleaning
- Supports both hygiene and visibility





The AT-1/AT-2 Applanation Tonometer features a seamless measuring prism that is easy to clean and suitable for repeated use. Designed to improve the accuracy of intraocular pressure measurement, it supports smooth clinical workflows.

## Measurement performance supporting routine clinical practice

Consistent measurement performance supports clinical decision-making.

Compliant with the Goldmann  
applanation tonometry method

Conforms to ISO standards  
(ISO 8612)



## Supports both hygiene and visibility

The prism is designed for hygiene control and repeated use, contributing to reliable measurements in routine clinical practice.

### Seamless structure designed for hygiene control

Seamless structure (this product)
Example of a segmented structure

Groove-free design integrating the contact and surrounding areas

The methacrylic resin prism provides high transparency and hardness for clear observation.

## Design focused on storage and operability

Compatible with a wide range of slit lamps and blends naturally with their appearance. When not in use, it can be moved aside while still attached to the slit lamp, allowing routine observation and enabling smooth transitions in daily examination workflows.

### [AT-1]

Compatible with a wide range of current TAKAGI slit lamps

### [AT-2]

Exclusively designed for 2ZL/4ZL slit lamps, providing a unified appearance when installed

| Specifications              |                                       |      |                                 |
|-----------------------------|---------------------------------------|------|---------------------------------|
| Model name                  | AT-1 Applanation Tonometer            |      | AT-2 Applanation Tonometer      |
| Attachable slit lamp models | 700GL, 700GL NSW, 700GL LPH, 4ZL, 2ZL | 30GL | 4ZL, 2ZL                        |
| Mounting tonometer base     | TB700                                 | TB30 | TB2                             |
| Measurement range           | 0-80 mmHg                             |      |                                 |
| Dimensions                  | 72 (W)×135 (D)×251 (H) mm             |      | 84.3 (W)×207.6 (D)×200.6 (H) mm |
| Weight                      | 0.44 kg                               |      | 0.67 kg                         |

\*The tonometer base is optional.

AT-1

AT-2

Learn more

To ensure safe and correct handling please read the user manual before using.

- Description and appearance as detailed in this brochure may be subject to change as improvements are made to products.
- Colors of product(s) in brochure may slightly differ from actual product(s) due to lighting situation when photographed or print conditions.

For The Americas, Asia-Pacific & Middle East  
**Takagi Seiko Co., Ltd.**  
 330-2 Iwafune, Nakano-shi, Nagano-ken, 383-8585, Japan  
 TEL : +81(0)269-22-6094 URL : <https://www.takagi-j.com>

For Europe & Africa  
**Takagi Ophthalmic Instruments Europe Ltd.**  
 Citylabs 1.0, Nelson Street, Manchester, M13 9NQ, UK  
 TEL : +44 (0)161 273 6330 URL : <https://www.takagieurope.com>

Takagi Seiko Co., Ltd.

1639 0120

Your local distributor

<https://www.takagi-global.com>

PB(AT-1/2)E Rev.0 Printed in Japan 2026.7 OV