

CR-750 CR-750S

Designed to support
clinical workflow


- 150 mm vertical stroke accommodates a wide range of patients
- Wide legroom and smooth operation
- Selectable seat colors to match clinic interiors



This motorized patient chair combines patient comfort with ease of operation for medical staff. Smooth height adjustment enables easy posture changes during examinations, while a design that blends with the clinic environment supports a cohesive clinical setting.

Patient comfort and efficient examinations

The 150 mm vertical stroke supports comfortable positioning for patients of various body types. In addition, the optional sliding mechanism allows patients to be positioned closer to examination equipment without moving the chair, helping to streamline examinations.



Staff-friendly operation

In addition to operation via the UP/DOWN switch on the back, the CR-750 can be operated using the standard foot controller, and the CR-750S can be operated from the ST-50 Smart Table. This allows for smooth and precise adjustments during use, reducing the workload for medical staff.

*The CR-750 cannot be connected to the ST-50.



Serving Your Vision

TAKAGI ophthalmic devices, proudly made in Japan since 1955.

Seat colors that complement clinic interiors

A choice of three seat colors* allows coordination with the clinic environment.

*The three seat colors are available only in the Americas, Asia, and the Middle East. Pink and green are not available in Europe and Africa.




Specifications

Model name	CR-750 Motorized Patient Chair CR-750S Motorized Patient Chair
Elevation stroke	150 mm
Load capacity	up to 150 kg
Input power	AC 100-230 V, 50-60 Hz
Power consumption	120 VA
Weight	26.5 kg (CR-750), 25 kg (CR-750S)

Dimensions

Dimensions (mm)



CR-750

Operable via foot controller



CR-750S

Controllable from the
ST-50 Smart Table



*The dimensions of CR-750 and CR-750S are identical.

*CR-750S does not include the foot controller.

Learn more



To ensure safe and correct handling please read the user manual before using.

- Description and appearance as detailed in this brochure may be subject to change as improvements are made to products.
- Colors of product(s) in brochure may slightly differ from actual product(s) due to lighting situation when photographed or print conditions.

For The Americas, Asia-Pacific & Middle East Takagi Seiko Co., Ltd.

330-2 Iwafune, Nakano-shi, Nagano-ken, 383-8585, Japan
TEL : +81(0)269-22-6094 URL : <https://www.takagi-j.com>

For Europe & Africa

Takagi Ophthalmic Instruments Europe Ltd.

Citylabs 1.0, Nelson Street, Manchester, M13 9NQ, UK
TEL : +44 (0)161 273 6330 URL : <https://www.takagi-europe.com>



Your local distributor